

企业简介

Corporation profile

Jiangxi Medwell medical products co., ltd located in the World Cultural Landscape World Geopark at the foot of Lushan Mountain, convenient transportation and elegant environment, 20 KM from the Lushan airport, 50 KM away from Nanchang Changbei airport, Waterway of the Changjiang, joint to chongqing, shanghai and oversea, which open up direct route from Jiujiang to Japan, Singapore, Philippines, Korea, Russian, Hong Kong, Taiwan and other countries. The international container liner open up Jiujiang to Shanghai through COSCO, SINOTRANS, APL, OOCL, CMA and PIL etc direct to all over the world.

THE company was established in 1991, through long-term exploration and unremitting efforts, introduction of advanced management concept, which has a good and healthy environment for production, advanced production equipment, complete facilities technology strong. Our main products are: with perfect quality products, good market reputation, the best cost-effective, our product exported to America, Europe, Asia and the Middle East.

The company has passed ISO13485:2003 certification and CE, to ensure the quality of products. Jiangxi Medwell medical products co., ltd. Adoption of the best hardware facilities and complete management system, a solid technological foundation and advanced modern enterprise management concepts together organically, always to provide customers with high-quality products and services.

Quality Guidelines: Adhere to the principle of prevention, ensure the quality meet customer requirements, positive assurance of product security and usability request in line with law and regulation.



MEDWELL MEDICAL PRODUCTS CO., LTD. 江西美宝利医用敷料有限公司

从质量到服务
我们将努力做的更好!

FROM QUALITY TO SERVICE
WE WILL DO BETTER

企业文化

Corporation culture

认真、诚实，树美宝利风貌。
交流，创新，建学习型企业。

Staidness, Honesty, Building up the mental
outlook of MEDWELL.

Intercourse, Innovation, Creating study type
enterprise.



企业文化

Corporation culture



MEETING ROOM 会议室

EXTERIOR SCENE OF THE PLANT 车间外景





企业精神

Corporation spirit

团结、忠诚、务实、学习。
Unity, Faithfulness, Practicableness, Study.

企业精神

Corporation spirit



质量方针

Corporation statement

坚持贯彻预防为主的方针，确保产品的质量满足顾客要求，保证产品的安全性和使用性能符合法规要求。
Combine prevention with elimination persistently
Ensure our product quality meeting customers' requirements and Guarantee the safety and usage of our products could tally with demands of rules and regulations.



质量方针

Corporation statement

WASHING HAND ROOM 洗手间

STERILE ROOM 消毒间



CHANGING CLOTHES ROOM 更衣室

WALKWAY 走道

生产车间
PRODUCTION WORKSHOP



WEAVING WORKSHOP

织布车间



PRODUCTION WORKSHOP
生产车间

生产车间
PRODUCTION WORKSHOP



AUTO MATIC FOLDING WORKSHOP



AUTO MATIC FOLDING WORKSHOP



生产车间
PRODUCTION WORKSHOP

生产车间
PRODUCTION WORKSHOP



KERLIK BANDAGE/SHRINKING BANDAGES
DRYING GAUZE ROLL WORKSHOP
曲缩烘干, 纱布卷车间



PRODUCTION WORKSHOP
生产车间

GAUZE BANDAGE 纱布绷带

- Gauze bandage is 100% made of bleaching absorbent cotton fabric with high absorbency and excellent permeability.
- It is ideal for wound packing and general wound care.
- Packing method: White Paper per roll, Blue kraft paper per dozen (12 rolls)



40'/30 x 20 White paper per roll, blue kraft per dozen (12)
40'/30 x 20 每卷白纸包 12卷蓝牛皮纸包装

Code	Size/Description	Amount/Unit	Units/Case
A40015Y	1"X5Y	12/Dozen	400/case
A40025Y	2"X5Y	12/Dozen	200/case
A40035Y	3"X5Y	12/Dozen	200/case
A40045Y	4"X5Y	12/Dozen	100/case
A40055Y	5"X5Y	12/Dozen	100/case
A40065Y	6"X5Y	12/Dozen	100/case
A40085Y	8"X5Y	12/Dozen	60/case
A40016Y	1"X6Y	12/Dozen	400/case
A40026Y	2"X6Y	12/Dozen	200/case
A40036Y	3"X6Y	12/Dozen	100/case
A40046Y	4"X6Y	12/Dozen	100/case
A40056Y	5"X6Y	12/Dozen	80/case
A40066Y	6"X6Y	12/Dozen	50/case
A40086Y	8"X6Y	12/Dozen	40/case
A40010M	1"X10m	12/Dozen	200/case
A400210M	2"X10m	12/Dozen	100/case
A400310M	3"X10m	12/Dozen	50/case
A400410M	4"X10m	12/Dozen	50/case
A400510M	5"X10m	12/Dozen	25/case
A400610M	6"X10m	12/Dozen	25/case
A400810M	8"X10m	12/Dozen	25/case



TRIANGULAR BANDAGE 三角绷带

Disposable Baby Receiving Blankets

Code	Size/Description	Amount/Unit	Units/Case
9500	32"x 40", 81.3cm x 101.6cm White with red and blue baby bottles	25/Pkgs	4/case

Triangular Bandages W/Safety Pins, Strong Muslin

Code	Size/Description	Amount/Unit	Units/Case
4910	40"x 40"x 56"	12/Bag	20/case
4915	36"x 36"x 51"	12/Bag	20/case
4916	36"x 36"x 51", Bleached	12/Bag	20/case

GAUZE ROLL 纱布卷

4-PLY Gauze Roll:

4-Ply gauze roll is 100% made of bleaching absorbent cotton fabric with high absorbency and soft touch without fluorescence. There are two types: With X-Ray threads and without X-Ray threads. Packing Method: White and Blue kraft paper.

2-PLY Gauze Roll:

2-Ply gauze roll is 100% made of bleaching absorbent cotton fabric with high absorbency and soft touch without fluorescence. There are two types: With X-Ray threads and without X-Ray threads. Packing Method: White and Blue kraft paper.

Open and Plain Gauze Roll:

Open and Plain gauze roll is 100% made of bleaching absorbent cotton fabric with high absorbency and soft touch without fluorescence. There are two types: With X-Ray and without X-Ray threads. Packing Method: Inner with white and blue kraft paper, outer polybag.

Fluff Bandage:

Open and Plain gauze roll is 100% made of bleaching absorbent cotton fabric with high absorbency and soft touch without fluorescence. There are two types: With X-Ray and without X-Ray threads. Packing Method: Inner with white and blue kraft paper, outer polybag.

Fluff Sponge:

Open and Plain gauze roll is 100% made of bleaching absorbent cotton fabric with high absorbency and soft touch without fluorescence. There are two types: With X-Ray and without X-Ray threads. Packing Method: Inner with white and blue kraft paper, outer polybag.



4-ply Gauze Roll Inner White Outer Blue Paper Pack
四折卷 内白外蓝纸包

40 S/30X20

Code	Size/Description	Amount/Unit	Units/Case
B4020091-4	90cmX91m	1roll	20/case
B4020100-4	90cmX100m	1roll	20/case
B4020050-4	90cmX50m	1roll	40/case
B4020005-4	90cmX5m	1roll	200/case

40 S/26X18

Code	Size/Description	Amount/Unit	Units/Case
B4017091-4	90cmX91m	1roll	20/case
B4017100-4	90cmX100m	1roll	20/case
B4017050-4	90cmX50m	1roll	40/case
B4017005-4	90cmX5m	1roll	200/case

40 S/19X15

Code	Size/Description	Amount/Unit	Units/Case
B4013091-4	90cmX91m	1roll	20/case
B4013100-4	90cmX100m	1roll	20/case
B4013050-4	90cmX50m	1roll	40/case
B4013005-4	90cmX5m	1roll	200/case

40 S/20X12

Code	Size/Description	Amount/Unit	Units/Case
B4012091-4	90cmX91m	1roll	20/case
B4012100-4	90cmX100m	1roll	20/case
B4012050-4	90cmX50m	1roll	40/case
B4012005-4	90cmX5m	1roll	200/case

40 S/12X8

Code	Size/Description	Amount/Unit	Units/Case
B4008091-4	90cmX91m	1roll	20/case
B4008100-4	90cmX100m	1roll	20/case
B4008050-4	90cmX50m	1roll	40/case
B4008005-4	90cmX5m	1roll	200/case


4-ply Gauze Roll Sewing Into X-ray , Inner white Outer Blue Paper Pack
4折卷 织入X光线 内白外蓝纸包装

Code	Size/Description	Amount/Unit	Units/Case
B4013091-4X	36*X100Y-4P	91m/roll	20/case
B4017091-4X	36*X100Y-4P	91m/roll	20/case
B4017100-4X	36*X100m-4P	100m/roll	20/case



2-ply Gauze Roll Inner white Outer Blue Paper Pack
对折纱布卷 内白外蓝纸包装

Code	Size/Description	Amount/Unit	Units/Case
C4020091-2	90cmX91m	1roll	20/case
C4020100-2	90cmX100m	1roll	20/case
C4017091-2	90cmX91m	1roll	20/case
C4017100-2	90cmX100m	1roll	20/case
C4013091-2	90cmX91m	1roll	20/case
C4013100-2	90cmX100m	1roll	20/case
C4012091-2	90cmX91m	1roll	20/case
C4012100-2	90cmX100m	1roll	20/case



平幅纱布卷：内白蓝纸，外聚酯膜包

Open And Plain Gauze Roll Inner With
And Blue Karft Paper Outer Polybag

Code	Size/Description	Amount/Unit	Units/Case
D40221000-A	40S/28X24 90cmX1000m	1roll	1roll
D40222000-A	40S/28X24 90cmX2000m	1roll	1roll
D40221000-B	40S/28X24 120cmX1000m	1roll	1roll
D40222000-B	40S/28X24 120cmX2000m	1roll	1roll
D40201000-A	40S/30X20 90cmX1000m	1roll	1roll
D40202000-A	40S/30X20 90cmX2000m	1roll	1roll
D40201000-B	40S/30X20 120cmX1000m	1roll	1roll
D40202000-B	40S/30X20 120cmX2000m	1roll	1roll
D40171000-A	40S/24X20 90cmX1000m	1roll	1roll
D40172000-A	40S/24X20 90cmX2000m	1roll	1roll
D40171000-B	40S/24X20 120cmX1000m	1roll	1roll
D40172000-B	40S/24X20 120cmX2000m	1roll	1roll
D40131000-A	40S/19X15 90cmX1000m	1roll	1roll
D40132000-A	40S/19X15 90cmX2000m	1roll	1roll
D40131000-B	40S/19X15 120cmX1000m	1roll	1roll
D40132000-B	40S/19X15 120cmX2000m	1roll	1roll



FLUFF BANDAGE FLUFF SPONGE

曲皱卷 曲缩片

40S/12X8 Fluff Bandage

Code	Size/Description	Amount/Unit	Units/Case
I4008-4.1Y	4.5"x4.1yds-6p	1roll	100/case
I4008-3.6Y	3.4"x3.6yds-6p	1roll	100/case

10S/12X8 Fluff Sponge, 32 Ply, Non-Sterile

Code	Size/Description	Amount/Unit	Units/Case
I4008-6.75"	6"x6.75"	100/bag	6/case
I4008-8.75"	7.75"x8.75"	100/bag	4/case
I4008-10.75"	10 1/2"x10 3/4"	50/bag	4/case

10S/12X8 Fluff Sponge, 32 Ply, Sterile

Code	Size/Description	Amount/Unit	Units/Case
I4008-6.75"	Medium, Sterile 2's, Soft Pouch	40/Pkgs	12/case
I4008-8.75"	Medium, Sterile 5's, Hard Pouch	5/Tray	120/case
I4008-10.75"	Medium, Sterile 10's, Hard Pouch	10/Tray	48/case

40S/28X24 RIBBON GAUZE BANDAGE , RIBBON GAUZE STRIP

内四折卷，内四折条

Code	Size/Description	Amount/Unit	Units/Case
J4022091-1/2"	1/2"x100yds-4p	1roll	10/case
J4022091-1"	1"x100yds-4p	1roll	10/case
J4022091-2"	2"x100yds-4p	1roll	10/case
J4022091-3"	3"x100yds-4p	1roll	10/case
J4022091-4"	4"x100yds-4p	1roll	10/case
J4022091-1.25cm	1.25cm x 100cm/4p	1roll	10/case
J4022091-1.5cm	1.5cm x 5.5cm/54p+x	1roll	10/case
J4022091-1.25cm	1.25cm x 50cm/4p	1roll	10/case
J4022091-1.25cm	1.25cm x 30cm/4p	1roll	10/case
J4022091-1.25cm	1.25cm x 10cm/4p	1roll	10/case
J4022091-1.25cm	1.25cm x 200cm/4p	1roll	10/case
J4022091-2cm	2cm x 75cm/4p	1roll	10/case



GAUZE SWABS 纱布片

- Gauze Sponge is 100% made of bleaching absorbent cotton fabric, with high absorbency and soft touch.
- With a wide range of specifications in size ply threads count and mesh.
- It is folded all by machine and perfect for absorbing blood and exudates.
- It is ideal for wound packing and general wound care.



GAUZE SWABS 纱布片

40S/19X15 Folded	40S/19X15 折边		
Code	Size/Description	Amount/Unit	Units/Case
E40132208-E	2"x2"-8p	100/pcs	50/case
E40132212-E	2"x2"-12p	100/pcs	50/case
E40133308-E	3"x3"-8p	100/pcs	40/case
E40133312-E	3"x3"-12p	100/pcs	40/case
E40134408-E	4"x4"-8p	100/pcs	40/case
E40134412-E	4"x4"-12p	100/pcs	20/case
E40134416-E	4"x4"-16p	100/pcs	20/case

40S/26X18 Folded	40S/26X18 折边		
Code	Size/Description	Amount/Unit	Units/Case
E40172208-E	2"x2"-8p	100/pcs	50/case
E40172212-E	2"x2"-12p	100/pcs	50/case
E40173308-E	3"x3"-8p	100/pcs	40/case
E40173312-E	3"x3"-12p	100/pcs	40/case
E40174408-E	4"x4"-8p	100/pcs	40/case
E40174412-E	4"x4"-12p	100/pcs	20/case
E40174416-E	4"x4"-16p	100/pcs	20/case

40S/20X12 Unfolded	40S/20X12 不折边		
Code	Size/Description	Amount/Unit	Units/Case
E40122208-C	2"x2"-8p	100/pcs	50/case
E40122212-C	2"x2"-12p	100/pcs	50/case
E40123308-C	3"x3"-8p	100/pcs	40/case
E40123312-C	3"x3"-12p	100/pcs	40/case
E40124408-C	4"x4"-8p	100/pcs	40/case
E40124412-C	4"x4"-12p	100/pcs	20/case
E40124416-C	4"x4"-16p	100/pcs	20/case



STERILE GAUZE SWABS 纱布片

40S/19X15 Folded, X-ray thread 40S/19X15 折边,X光线

Code	Size/Description	Amount/Unit	Units/Case
E40132208-EX	2"x2"-8p	100/pcs	50/case
E40132212-EX	2"x2"-12p	100/pcs	50/case
E40133308-EX	3"x3"-8p	100/pcs	40/case
E40133312-EX	3"x3"-12p	100/pcs	40/case
E40134408-EX	4"x4"-8p	100/pcs	40/case
E40134412-EX	4"x4"-12p	100/pcs	20/case
E40134416-EX	4"x4"-16p	100/pcs	20/case

40S/26X18 Folded, X-ray thread 40S/26X18 折边,X光线

Code	Size/Description	Amount/Unit	Units/Case
E40172208-EX	2"x2"-8p	100/pcs	50/case
E40172212-EX	2"x2"-12p	100/pcs	50/case
E40173308-EX	3"x3"-8p	100/pcs	40/case
E40173312-EX	3"x3"-12p	100/pcs	40/case
E40174408-EX	4"x4"-8p	100/pcs	40/case
E40174412-EX	4"x4"-12p	100/pcs	20/case
E40174416-EX	4"x4"-16p	100/pcs	20/case

40S/20X12 Unfolded, X-ray thread 40S/20X12 不折边,X光线

Code	Size/Description	Amount/Unit	Units/Case
E40122208-CX	2"x2"-8p	100/pcs	50/case
E40122212-CX	2"x2"-12p	100/pcs	50/case
E40123308-CX	3"x3"-8p	100/pcs	40/case
E40123312-CX	3"x3"-12p	100/pcs	40/case
E40124408-CX	4"x4"-8p	100/pcs	40/case
E40124412-CX	4"x4"-12p	100/pcs	20/case
E40124416-CX	4"x4"-16p	100/pcs	20/case



STERILE MOLD PACKING 吸塑包装

40S/20X12 Unfolded, 5'S Mold packing 40S/20X12不折边 5片 吸塑包装

Code	Size/Description	Amount/Unit	Units/Case
G40122208-C5	2"x2"-8p	5/pkg 20/box	20/case
G40122212-C5	2"x2"-12p	5/pkg 20/box	20/case
G40123308-C5	3"x3"-8p	5/pkg 20/box	15/case
G40123312-C5	3"x3"-12p	5/pkg 20/box	15/case
G40124408-C5	4"x4"-8p	5/pkg 20/box	10/case
G40124412-C5	4"x4"-12p	5/pkg 20/box	10/case
G40124416-C5	4"x4"-16p	5/pkg 20/box	20/case

40S/20X12 Folded, 5'S Mold packing 40S/20X12折边 5片 吸塑包装

Code	Size/Description	Amount/Unit	Units/Case
F40172208-E5	2"x2"-8p	5/pkg 20/box	20/case
F40172212-E5	2"x2"-12p	5/pkg 20/box	20/case
F40173308-E5	3"x3"-8p	5/pkg 20/box	15/case
F40173312-E5	3"x3"-12p	5/pkg 20/box	15/case
F40174408-E5	4"x4"-8p	5/pkg 20/box	10/case
F40174412-E5	4"x4"-12p	5/pkg 20/box	10/case
F40174416-E5	4"x4"-16p	5/pkg 20/box	20/case



STERILE MOLD PACKING
消毒吸垫包装

40S/19X15 Folded, 1'S sterile packing 40S/19X15 折边 1片 消毒袋包装

Code	Size/Description	Amount/Unit	Units/Case
F40132208-E1	2"x2"-8p	1/pouch 100/box	24/case
F40132212-E1	2"x2"-12p	1/pouch 100/box	24/case
F40133308-E1	3"x3"-8p	1/pouch 100/box	24/case
F40133312-E1	3"x3"-12p	1/pouch 100/box	24/case
F40134408-E1	4"x4"-8p	1/pouch 100/box	12/case
F40134412-E1	4"x4"-12p	1/pouch 100/box	12/case

40S/26X18 Folded, 1'S sterile packing 40S/26X18 折边 1片 消毒袋包装

Code	Size/Description	Amount/Unit	Units/Case
F40172208-E1	2"x2"-8p	1/pouch 100/box	24/case
F40172212-E1	2"x2"-12p	1/pouch 100/box	24/case
F40173308-E1	3"x3"-8p	1/pouch 100/box	24/case
F40173312-E1	3"x3"-12p	1/pouch 100/box	24/case
F40174408-E1	4"x4"-8p	1/pouch 100/box	12/case
F40174412-E1	4"x4"-12p	1/pouch 100/box	12/case
F40174416-E1	4"x4"-16p	1/pouch 100/box	12/case

40S/26X18 Unfolded, 2'S sterile packing

40S/26X18 不折边 2片 消毒袋包装

Code	Size/Description	Amount/Unit	Units/Case
F40122208-C2	2"x2"-8p	2/pouch 50/box	24/case
F40122212-C2	2"x2"-12p	2/pouch 50/box	24/case
F40123308-C2	3"x3"-8p	2/pouch 50/box	24/case
F40123312-C2	3"x3"-12p	2/pouch 50/box	24/case
F40124408-C2	4"x4"-8p	2/pouch 50/box	12/case
F40124412-C2	4"x4"-12p	2/pouch 50/box	12/case
F40124416-C2	4"x4"-16p	2/pouch 50/box	12/case

GAUZE SPONGES U-TYPE
开口片 (纱布, 无纺布)

Specification

7.5cm x 7.5cm	8ply, 12ply, 16ply
10cm x 10cm	8ply, 12ply, 16ply
10cm x 20cm	8ply, 12ply, 16ply
10cm x 10cm	24ply, 32ply, 64ply



NON-WOVEN SPONGES
无纺布片

Specification

5cm x 5cm	8ply, 12ply, 16ply
5cm x 7.5cm	8ply, 12ply, 16ply
7.5cm x 7.5cm	8ply, 12ply, 16ply



Non-sterile packing Sterile Packing

*spnlace non-woven, medical grade, high absorbency and soft
*70% visconse, 30% polyester
*weight: 30g/cm2, 35g/cm2, 40g/cm2 or 50g/cm2

❁ LAP SPONGE 腹部垫

- Lap Sponge is 100% made of bleaching absorbent cotton fabric, with high absorbency.
- It is offered to O.R. team for absorbing blood and exudates or for cleaning. It can prevent.
- The wound from being effect.
- The sterile methods has two types: Steam and Gamma Ray.
- Our Lap Sponge is sewn and cross-stitched.



➔ LAP SPONGE 腹部垫

40S/24 x 20 Washed Blue loop, X-ray chip, 5'S
40S/24 x 20 蓝带、锁片、水洗 5片/扎

Code	Size/Description	Amount/Unit	Units/Case
H40174545-2	45x45cm-2p	5/pkg	320/case
H40174545-4	45x45cm-4p	5/pkg	160/case
H40174545-6	45x45cm-6p	5/pkg	120/case
H40174545-8	45x45cm-8p	5/pkg	100/case
H40174040-2	40x40cm-2p	5/pkg	560/case
H40174040-4	40x40cm-4p	5/pkg	280/case
H40174040-6	40x40cm-6p	5/pkg	200/case
H40174040-8	40x40cm-8p	5/pkg	160/case

40S/24 x 20 Blue loop, X-ray chip, 5'S packing
by CSR paper Mold packing
40S/24 x 20 蓝带、锁片、5片/CSR包 吸塑包装

Code	Size/Description	Amount/Unit	Units/Case
H40174545-4A	45x45cm-4p	5/pkg	20/case
H40174040-4A	40x40cm-4p	5/pkg	20/case
H40173030-4A	30x30cm-4p	5/pkg	20/case

40S/24 x 20 Blue loop, X-ray chip, 5'S packing
by wrinkle paper, Mold packing
40S/24 x 20 蓝带、锁片、5片/皱纸包 吸塑包装

Code	Size/Description	Amount/Unit	Units/Case
H40174545-4B	45x45cm-4p	5/pkg	20/case
H40174040-4B	40x40cm-4p	5/pkg	20/case
H40173030-4B	30x30cm-4p	5/pkg	20/case



→ LAP SPONGE 腹巾

40S/19 x 15 Unwashed Blue loop, X-ray chip, 5'S
40S/19 x 15 蓝带、钢片、不水洗 5片/扎

Code	Size/Description	Amount/Unit	Units/Case
H40134545-2	45x45cm-2p	5/pkg	320/case
H40134545-4	45x45cm-4p	5/pkg	160/case
H40134545-6	45x45cm-6p	5/pkg	120/case
H40134545-8	45x45cm-8p	5/pkg	100/case
H40134040-2	40x40cm-2p	5/pkg	560/case
H40134040-4	40x40cm-4p	5/pkg	280/case
H40134040-6	40x40cm-6p	5/pkg	200/case
H40134040-8	40x40cm-8p	5/pkg	160/case

✿ GAUZE BALL 纱布球

- Gauze ball is 100% made of bleaching absorbent cotton fabric with high absorbency.
- There are two types: With X-Ray threads and without X-Ray threads.
- It is mainly used to absorbing blood and exudates during the operation and cleaning the wounds.



GAUZE BALL 纱布球

40S/30 x 20 Round Gauze Ball with X-ray and non X-ray
40S/30 x 20 X光线与不带X光线纱布球

Code	Size/Description	Amount/Unit	Units/Case
K40201010	10x10cm	250/polybag	32/case
K40201515	15x15cm	250/polybag	24/case
K40202020	20x20cm	250/polybag	16/case
K40203030	30x30cm	250/polybag	8/case

40S/26 x 18 Round Gauze Ball with X-ray and non X-ray
40S/26 x 18 X光线与不带X光线纱布球

Code	Size/Description	Amount/Unit	Units/Case
K40171010	10x10cm	250/polybag	32/case
K40171515	15x15cm	250/polybag	24/case
K40172020	20x20cm	250/polybag	16/case
K40173030	30x30cm	250/polybag	8/case
K40174040	40x40cm	250/polybag	8/case

40S/26 x 18 Round Gauze, Mold packing
40S/26 x 18 纱布球, 模塑包装

Code	Size/Description	Amount/Unit	Units/Case
K40172020-2	20x20cm	2/pkg 15/box	60/case
K40172020-3	20x20cm	2/pkg 15/box	45/case
K40172020-5	20x20cm	2/pkg 15/box	45/case

20 SINGLE COTTON TIPPED APPLICATORS 棉签

- The CTA'S are packaged in convenient, peel-down pouches and are available in sterile and non-sterile, wooden or plastic shafts.
- One of our latest product offerings is the 100% cotton, single tipped applicators in a variety of sizes and put-ups.